

Mindray DP-5

Hersteller: Mindray

Art.Nr.: MI_9950B-PA00012

Produktbeschreibung

15-inch LCD Monitor Three Transducer Connectors iBeam™ (Spatial Compounding Imaging) iClear™ (Speckle Suppression Imaging) iTouch™ (Auto Image Optimization) PSH™ (Phase Shift Harmonic Imaging) Trapezoid Imaging/ExFOV Imaging Modes: B/2B/4B/M/B+M 8-segment TGC 320GB Hard Disk & iStation™ Patient Information Management Shared Service (includes Obstetrics, Pediatrics, Gynecology, Abdomen, Urology, Small Parts, Vascular, Cardiology, Nerve and Emergency Medicine packages?)

Technische Daten

3D Darstellung	—	Akkubetrieb	○
Anzahl aktiver Sondenports	3	Anzahl USB Ports	2
Auto Nackenfaltenmessung	○	Automatische Kurvenvermessung	—
Automatische Livedarstellung	—	Autooptimierung	●
Bluetooth	○	Color Flow Modus	—
Compound Modus	●	CW Doppler	—
DIOCM mit Worklist	○	Duplexmodus	—
DVD-Laufwerk	●	Echtzeit 4D Darstellung	—
Elastographie	—	Erweiterter Rauschfilter	●
Farbgerät	—	Festplatte	●
FullScreen Bilddarstellung	—	Gewebespezifischer Filter	●
Harmonic Imaging	—	Kalk. Geburtshilfe OB	●
Kalk. Gyn	●	Kalk. Renal	●
Kalk. Uro	●	Kalkulationsprogramm Cardio	—
Kalkulationsprogramm Vaskular	—	Konsolengerät	●
Monitorgröße	15"	Netzwerkspeicherung	●
Panoramaview	—	Parameter nach Freeze verstellbar	—
Patientendatenbank	●	PDI	—
PW Doppler	—	Rauschfilter	●
Touchscreen	—	Tragbar	—
Trapezoidmodus	●	Triplex Modus	—
TVI	—	Videoprinterslot	●
Videoschnittstellen	PAL, S-VHS, DVI	Wlan	○
Zoomfunktion	●		

verfügbare Optionen

- ▶ **Mindray iPower (built-in batter, supports working time for 1.5 hour) (Art.Nr. MI_iPower_DP5)**
iPower (built-in batter, supports working time for 1.5 hour)

kompatible Produkte

- ▶ **Mindray 35C20EA Elektronischer Mikrokonvex-Array-Schallkopf (Art.Nr. MI_35C20EA)**
Mikrokonvexsonde für 2,0/3,5/6,0 MHz, Abdominal, OB/Gyn, Vet,
- ▶ **Mindray 35C50EA Elektronischer Konvex-Array Schallkopf (Art.Nr. MI_35C50EA)**
Abdomenschallkopf 2.5/3.5/5.0 MHz für DP-6600, DP-6900, DP-8500, DP-50, DP-5, Abdominal, OB/Gyn
- ▶ **Mindray 65C15EA Elektronischer Mikrokonvex-Array-Schallkopf (Art.Nr. MI_65C15EA)**
Elektronischer Mikrokonvex-Array-Schallkopf, 5,0/6,5/8,0 MHz, Abdomen, Kleinkinder
- ▶ **Mindray 65EB10EA elektr. biplanare Sonde (Art.Nr. MI_65EB10EA)**
Schallsonde biplanar für die Rektaluntersuchung, 5.0 - 8.5 MHz. Leider noch kein Bild verfügbar.
- ▶ **Mindray 65EC10EA Endokavitär-Schallkopf (Art.Nr. MI_65EC10EA)**
Elektronischer Endokavität-Schallkopf (5,0/6,5/8,0 MHz), Transvaginal, Transrectal
- ▶ **Mindray 75L38EA Elektronischer Linear-Array-Schallkopf (Art.Nr. 75L38EA)**
Linearsonde 38mm, 5,0 7,5 10,0 MHz, Vascular, Small Parts, MSK
- ▶ **Mindray 75L53EA Elektronischer Linear-Array-Schallkopf (Art.Nr. MI_75L53EA)**
Linearsonde für die Geräte DP-50 und DP-5, 5,0/7,5/10,0 MHz, Vascular, Small Parts, MSK
- ▶ **Mindray Auto IMT Package (Automatic measurement for Intima-Media Thickness) (Art.Nr. MI_110-001962-00)**
Auto IMT Package (Automatic measurement for Intima-Media Thickness)
- ▶ **Mindray DICOM Basic (Art.Nr. MI_110-001963-00)**
DICOM Basic
- ▶ **Mindray DICOM Worklist (Art.Nr. MI_110-001964-00)**
DICOM Worklist
- ▶ **Mindray MITSUBISHI P93W (analog B/W Printer) (Art.Nr. MI_0000-10-11057)**
MITSUBISHI P93W (analog B/W Printer)
- ▶ **Mindray SONY UP-X898MD (Digital/Analog printer with EU power cord, include BNC connector) (Art.Nr. MI_115-032024-00)**
SONY UP-X898MD (Digital/Analog printer with EU power cord, include BNC connector)
- ▶ **Mindray Water-resistant foot switch with three hot keys (USB port) (Art.Nr. MI_115-005252-00)**
Water-resistant foot switch with three hot keys (USB port)
- ▶ **Mindray Water-resistant foot switch with two hot keys (USB port) (Art.Nr. MI_0020-30-90810)**
Water-resistant foot switch with two hot keys (USB port)